

Product datasheet for **GA106652**

Human PSMD6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD6 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 6
Symbol:	PSMD6
UniProt ID:	Q15008
Synonyms:	p42A; p44S10; Rpn7; S10; SGA-113M
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)