

Product datasheet for GA106672

Human LRRC37A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC37A CRISPRa kit – CRISPR gene activation of human leucine rich repeat containing 37A
Symbol:	LRRC37A
Components:	GA106672G1, LRRC37A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106672G2, LRRC37A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106672G3, LRRC37A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A6NMS7</u>
Synonyms:	LRRC37
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

