

Product datasheet for **GA106677**

Human LPPR4 (PLPPR4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLPPR4 CRISPRa kit - CRISPR gene activation of human phospholipid phosphatase related 4
Symbol:	PLPPR4
Components:	GA106677G1, LPPR4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106677G2, LPPR4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106677G3, LPPR4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z2D5</u>
Synonyms:	LPR4; PHP1; PRG-1; PRG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)