

Product datasheet for **GA106680**

Human **TELO2** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TELO2 CRISPRa kit - CRISPR gene activation of human telomere maintenance 2
Symbol:	TELO2
Components:	GA106680G1, TELO2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106680G2, TELO2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106680G3, TELO2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y4R8</u>
Synonyms:	CLK2; TEL2; YHFS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)