

Product datasheet for **GA106693**

Human SPG48 (AP5Z1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AP5Z1 CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 5 subunit zeta 1
Symbol:	AP5Z1
Components:	GA106693G1, SPG48 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106693G2, SPG48 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106693G3, SPG48 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q43299</u>
Synonyms:	KIAA0415; SPG48; zeta
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)