

Product datasheet for **GA106698**

Human RICH2 (ARHGAP44) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARHGAP44 CRISPRa kit - CRISPR gene activation of human Rho GTPase activating protein 44
Symbol:	ARHGAP44
Components:	GA106698G1, RICH2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106698G2, RICH2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106698G3, RICH2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q17R89</u>
Synonyms:	NPC-A-10; RICH2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)