

Product datasheet for **GA106703**

Human CNAP1 (NCAPD2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NCAPD2 CRISPRa kit - CRISPR gene activation of human non-SMC condensin I complex subunit D2
Symbol:	NCAPD2
Components:	GA106703G1, CNAP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106703G2, CNAP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106703G3, CNAP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q15021
Synonyms:	CAP-D2; CNAP1; hCAP-D2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)