

Product datasheet for **GA106716**

Human KIAA0020 (PUM3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PUM3 CRISPRa kit - CRISPR gene activation of human pumilio RNA binding family member 3
Symbol:	PUM3
UniProt ID:	Q15397
Synonyms:	HBV X-transactivated gene 5 protein; HLA-HA8; KIAA0020 protein; MGC8749; minor histocompatibility antigen HA-8; PEN; penguin homolog; protein 5 transactivated by hepatitis B virus X antigen (HBxAg); PUF6; XTP5; XTP5, HLA-HA8, MGC8749, RP11-526D20.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)