

## Product datasheet for **GA106717**

### Human GPR105 (P2RY14) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2RY14 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2Y14
Symbol:	P2RY14
UniProt ID:	<a href="#">Q15391</a>
Synonyms:	BPR105; GPR105; P2Y14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)