

Product datasheet for **GA106733**

Human HS3ST4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HS3ST4 CRISPRa kit - CRISPR gene activation of human heparan sulfate-glucosamine 3-sulfotransferase 4
Symbol:	HS3ST4
Components:	GA106733G1, HS3ST4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106733G2, HS3ST4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106733G3, HS3ST4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y661</u>
Synonyms:	3-OST-4; 3OST4; 30ST4; h3-OST-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)