

Product datasheet for **GA106735**

Human HS3ST3A1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HS3ST3A1 CRISPRa kit - CRISPR gene activation of human heparan sulfate-glucosamine 3-sulfotransferase 3A1
Symbol:	HS3ST3A1
UniProt ID:	Q9Y663
Synonyms:	3-OST-3A; 3OST3A1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)