

Product datasheet for **GA106739**

Human USP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	USP3 CRISPRa kit - CRISPR gene activation of human ubiquitin specific peptidase 3
Symbol:	USP3
UniProt ID:	Q9Y6I4
Synonyms:	SIH003; UBP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)