

## Product datasheet for **GA106744**

### Human TL1A (TNFSF15) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNFSF15 CRISPRa kit - CRISPR gene activation of human TNF superfamily member 15
Symbol:	TNFSF15
UniProt ID:	<a href="#">O95150</a>
Synonyms:	TL1; TL1A; TNLG1B; VEGI; VEGI192A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)