

Product datasheet for **GA106755**

Human DOP1B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DOP1B CRISPRa kit - CRISPR gene activation of human DOP1 leucine zipper like protein B
Symbol:	DOP1B
UniProt ID:	Q9Y3R5
Synonyms:	21orf5; C21orf5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)