

Product datasheet for **GA106762**

Human PP4R1 (PPP4R1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP4R1 CRISPRa kit - CRISPR gene activation of human protein phosphatase 4 regulatory subunit 1
Symbol:	PPP4R1
UniProt ID:	Q8TF05
Synonyms:	MEG1; PP4(Rmeg); PP4R1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)