

## Product datasheet for **GA106777**

### Human **KCNE3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNE3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily E regulatory subunit 3
Symbol:	KCNE3
UniProt ID:	<a href="#">Q9Y6H6</a>
Synonyms:	HOKPP; HYPP; MiRP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)