

Product datasheet for **GA106787**

Human LNK (SH2B3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH2B3 CRISPRa kit - CRISPR gene activation of human SH2B adaptor protein 3
Symbol:	SH2B3
UniProt ID:	Q9UQQ2
Synonyms:	IDDM20; LNK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)