

## Product datasheet for **GA106789**

### Human HCN4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HCN4 CRISPRa kit - CRISPR gene activation of human hyperpolarization activated cyclic nucleotide gated potassium channel 4
Symbol:	HCN4
UniProt ID:	<a href="#">Q9Y3Q4</a>
Synonyms:	SSS2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)