

Product datasheet for **GA106794**

Human PIGK activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIGK CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class K
Symbol:	PIGK
UniProt ID:	Q92643
Synonyms:	GPI8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)