

Product datasheet for **GA106797**

Human PARP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PARP3 CRISPRa kit - CRISPR gene activation of human poly(ADP-ribose) polymerase family member 3
Symbol:	PARP3
UniProt ID:	<u>Q9Y6F1</u>
Synonyms:	ADPRT3; ADPRTL2; ADPRTL3; ARTD3; IRT1; PADPRT-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)