

Product datasheet for **GA106802**

Human SH2D3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH2D3A CRISPRa kit - CRISPR gene activation of human SH2 domain containing 3A
Symbol:	SH2D3A
Components:	GA106802G1, SH2D3A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106802G2, SH2D3A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106802G3, SH2D3A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BRG2</u>
Synonyms:	NSP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)