

Product datasheet for **GA106808**

Human SMC4 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | SMC4 CRISPRa kit - CRISPR gene activation of human structural maintenance of chromosomes 4 |
| Symbol: | SMC4 |
| UniProt ID: | Q9NTJ3 |
| Synonyms: | CAPC; CAPC, SMC4L1, hCAP-C; chromosome-associated polypeptide C; hCAP-C; SMC4 structural maintenance of chromosomes 4-like 1; SMC4L1; structural maintenance of chromosomes 4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)