

Product datasheet for **GA106826**

Human DPP3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | DPP3 CRISPRa kit - CRISPR gene activation of human dipeptidyl peptidase 3 |
| Symbol: | DPP3 |
| UniProt ID: | <u>Q9NY33</u> |
| Synonyms: | DPPIII |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)