

Product datasheet for **GA106830**

Human PHEMX (TSPAN32) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPAN32 CRISPRa kit - CRISPR gene activation of human tetraspanin 32
Symbol:	TSPAN32
UniProt ID:	<u>Q96QS1</u>
Synonyms:	MGC22455; PHEMX; PHMX; TSSC6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)