

Product datasheet for **GA106833**

Human PDCD7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDCD7 CRISPRa kit – CRISPR gene activation of human programmed cell death 7
Symbol:	PDCD7
Components:	GA106833G1, PDCD7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106833G2, PDCD7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106833G3, PDCD7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N8D1</u>
Synonyms:	ES18; HES18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

