

## Product datasheet for **GA106840**

### Human **KCNK7** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNK7 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 7
Symbol:	KCNK7
UniProt ID:	<a href="#">Q9Y2U2</a>
Synonyms:	K2p7.1; TWIK3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)