

Product datasheet for **GA106848**

Human TSPAN5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPAN5 CRISPRa kit - CRISPR gene activation of human tetraspanin 5
Symbol:	TSPAN5
UniProt ID:	P62079
Synonyms:	NET-4; NET4; TM4SF9; TSPAN-5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)