


## Product datasheet for **GA106849**

### Human TSPAN3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TSPAN3 CRISPRa kit - CRISPR gene activation of human tetraspanin 3
Symbol:	TSPAN3
UniProt ID:	<a href="#">O60637</a>
Synonyms:	TM4-A; TM4SF8; TSPAN-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)