

Product datasheet for **GA106862**

Human HIPK3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIPK3 CRISPRa kit - CRISPR gene activation of human homeodomain interacting protein kinase 3
Symbol:	HIPK3
UniProt ID:	<u>Q9H422</u>
Synonyms:	DYRK6; FIST3; PKY; YAK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)