

## Product datasheet for **GA106863**

### Human FEM 1B (FEM1B) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | FEM1B CRISPRa kit - CRISPR gene activation of human fem-1 homolog B                                                                                                                                                                                                                                           |
| Symbol:       | FEM1B                                                                                                                                                                                                                                                                                                         |
| Components:   | GA106863G1, FEM 1B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA106863G2, FEM 1B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA106863G3, FEM 1B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9UK73</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | F1A-ALPHA; F1AA; FEM1-beta                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)