

Product datasheet for **GA106867**

Human ARL4C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARL4C CRISPRa kit - CRISPR gene activation of human ADP ribosylation factor like GTPase 4C
Symbol:	ARL4C
Components:	GA106867G1, ARL4C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106867G2, ARL4C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106867G3, ARL4C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P56559</u>
Synonyms:	ARL7; LAK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)