

## Product datasheet for **GA106870**

### Human Dynein light chain (DNAL4) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	DNAL4 CRISPRa kit - CRISPR gene activation of human dynein axonemal light chain 4
Symbol:	DNAL4
UniProt ID:	<a href="#">O96015</a>
Synonyms:	MRMV3; PIG27
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)