

Product datasheet for **GA106875**

Human Hsp75 (TRAP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRAP1 CRISPRa kit - CRISPR gene activation of human TNF receptor associated protein 1
Symbol:	TRAP1
UniProt ID:	Q12931
Synonyms:	HSP 75; HSP75; HSP90L; TRAP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)