

Product datasheet for **GA106879**

Human Elastase 3A (CELA3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CELA3A CRISPRa kit - CRISPR gene activation of human chymotrypsin like elastase 3A
Symbol:	CELA3A
UniProt ID:	<u>P09093</u>
Synonyms:	ELA3; ELA3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)