

Product datasheet for **GA106884**

Human AKAP9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AKAP9 CRISPRa kit - CRISPR gene activation of human A-kinase anchoring protein 9
Symbol:	AKAP9
Components:	GA106884G1, AKAP9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106884G2, AKAP9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106884G3, AKAP9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99996
Synonyms:	AKAP-9; AKAP350; AKAP450; CG-NAP; HYPERION; LQT11; MU-RMS-40.16A; PPP1R45; PRKA9; YOTIAO
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)