

## Product datasheet for GA106901

### Human P2RY5 (LPAR6) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LPAR6 CRISPRa kit - CRISPR gene activation of human lysophosphatidic acid receptor 6
Symbol:	LPAR6
UniProt ID:	<a href="#">P43657</a>
Synonyms:	ARWH1; HYPT8; LAH3; P2RY5; P2Y5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

