

Product datasheet for **GA106927**

Human PSME3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSME3 CRISPRa kit - CRISPR gene activation of human proteasome activator subunit 3
Symbol:	PSME3
UniProt ID:	P61289
Synonyms:	HEL-S-283; Ki; PA28-gamma; PA28G; PA28gamma; REG-GAMMA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)