

## Product datasheet for **GA106964**

### Human hnRNP R (HNRNPR) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	HNRNPR CRISPRa kit - CRISPR gene activation of human heterogeneous nuclear ribonucleoprotein R
Symbol:	HNRNPR
Components:	GA106964G1, hnRNP R gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106964G2, hnRNP R gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106964G3, hnRNP R gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q43390</a></u>
Synonyms:	hnRNP-R; HNRPR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)