

Product datasheet for **GA106965**

Human **SLC35B1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35B1 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member B1
Symbol:	SLC35B1
UniProt ID:	P78383
Synonyms:	UGTREL1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)