

Product datasheet for **GA106967**

Human AP3S2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AP3S2 CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 3 subunit sigma 2
Symbol:	AP3S2
UniProt ID:	P59780
Synonyms:	AP3S3; sigma3b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)