

Product datasheet for **GA106983**

Human CNKSR1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CNKSR1 CRISPRa kit - CRISPR gene activation of human connector enhancer of kinase suppressor of Ras 1 |
| Symbol: | CNKSR1 |
| UniProt ID: | Q969H4 |
| Synonyms: | CNK; CNK1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)