

## Product datasheet for **GA107047**

### Human PCGF3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PCGF3 CRISPRa kit - CRISPR gene activation of human polycomb group ring finger 3
Symbol:	PCGF3
UniProt ID:	<a href="#">Q3KNV8</a>
Synonyms:	DONG1; RNF3; RNF3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)