

## Product datasheet for **GA107091**

### Human HEC1 (NDC80) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | NDC80 CRISPRa kit - CRISPR gene activation of human NDC80 kinetochore complex component                                                                                                                                                                                                                 |
| Symbol:       | NDC80                                                                                                                                                                                                                                                                                                   |
| Components:   | GA107091G1, HEC1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA107091G2, HEC1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA107091G3, HEC1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q14777</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | HEC; HEC1; HsHec1; hsNDC80; KNTC2; TID3                                                                                                                                                                                                                                                                 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)