

Product datasheet for **GA107096**

Human IFITM3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IFITM3 CRISPRa kit - CRISPR gene activation of human interferon induced transmembrane protein 3
Symbol:	IFITM3
UniProt ID:	<u>Q01628</u>
Synonyms:	1-8U; DSPA2b; IP15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)