

Product datasheet for **GA107108**

Human ARIH2 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | ARIH2 CRISPRa kit - CRISPR gene activation of human ariadne RBR E3 ubiquitin protein ligase 2 |
| Symbol: | ARIH2 |
| UniProt ID: | <u>O95376</u> |
| Synonyms: | ARI2; TRIAD1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)