

Product datasheet for **GA107155**

Human PCID1 (EIF3M) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | EIF3M CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit M |
| Symbol: | EIF3M |
| UniProt ID: | Q7L2H7 |
| Synonyms: | B5; GA17; hfl-B5; PCID1; TANGO7 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)