

Product datasheet for **GA107158**

Human SEC23B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEC23B CRISPRa kit - CRISPR gene activation of human SEC23 homolog B, coat complex II component
Symbol:	SEC23B
UniProt ID:	Q15437
Synonyms:	CDA-II; CDAII; CDAN2; HEMPAS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)