

Product datasheet for **GA107159**

Human SEC23 (SEC23A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEC23A CRISPRa kit - CRISPR gene activation of human Sec23 homolog A, coat complex II component
Symbol:	SEC23A
UniProt ID:	Q15436
Synonyms:	CLSD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)