

Product datasheet for **GA107167**

Human hnRNP Q (SYNCRIP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYNCRIP CRISPRa kit - CRISPR gene activation of human synaptotagmin binding cytoplasmic RNA interacting protein
Symbol:	SYNCRIP
UniProt ID:	<u>O60506</u>
Synonyms:	GRY-RBP; GRYRBP; hnRNP-Q; HNRNPQ; HNRNPQ1; NSAP1; PP68
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)