

Product datasheet for **GA107213**

Human TIP30 (HTATIP2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HTATIP2 CRISPRa kit - CRISPR gene activation of human HIV-1 Tat interactive protein 2
Symbol:	HTATIP2
UniProt ID:	<u>Q9BUP3</u>
Synonyms:	CC3; SDR44U1; TIP30
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)